

Offer patients the opportunity for better eye health and vision with proactive myopia management.






Lenses for Orthokeratology & Myopia Control

CORRECTING THE VISION OF YOUR PATIENT WHILE THEY SLEEPING

There is no need of an expensive trial lens set Effectiveness with the first lens





MOLD THE WAY YOUR PATIENTS SEE

Give your patients crisp, clear vision without the typical options of glasses, daytime contact lenses, or LASIK surgery.

Grow your practice while thrilling your patients with a proven vision correction option that can improve their lives.

CLINICALLY PROVEN EMPIRICAL FIT OUTSTANDING RESULTS


Emerald contact lenses are FDA approved and provide easy vision correction by gently reshaping the cornea at night while the patient sleeps.

They are safe, reversible and have been successfully adapted to over a million eyes worldwide with excellent results. Empirically adapted. Clinically proven. With surprising results.

SIMPLIFIED ADVANCED CORNEAL MOLDING

- Proven success rate of 85% in the first adaptation.
- \bullet Proven 95.5% success rate with up to 1 change.
- Clinically proven to reduce the progression of myopia in children.
- Simple adaptation only with refraction, keratometry and HVID.
- Empirically adapted no trial lens required.
- Clinically proven to suppress axial length elongation in childhood myopia.





UNIQUE MULTI-ZONE DESIGN

RADIO OF THE REAR OPTICAL ZONE	Designed to reshape the cornea to the required correction. Unlike a normal GP lens.
REVERSE CURVE	Determined by a clinically proven formula that aligns the base curve with the peripheral cornea and controls the apical clearance of the radius of the posterior optic zone.
ALIGNMENT AREA	A series of two or more separate curves, this zone determines the characteristics of the lens fit. It is designed to align closely with the cornea and provide proper centering for maximum performance.
PERIPHERAL CURVE	Peripheral curve provides clearance at the edge for proper tear exchange and easy lens removal.

Make a Difference in the Myopia Epidemic

Myopia is more than a refractive condition. It's a global epidemic that's having a serious impact on today's children around the world.

The prevalence of myopia has nearly doubled since the 1970s, ¹ and it is expected to become the leading cause of permanent blindness.

5 billion people or

of the world's population will be myopic by 2050



Myopia significantly increases the risk of more serious ocular health issues, including:

Glaucoma

14 4 X


for myopia >8.00D

Cataract

3 3 X

for myopia >6.00D

7.8 X





Advancing the **Geometry of Sight**

Euclid, the Greek mathematician considered the father of geometry, inspired the creation of the Euclid Emerald myopia management treatment.

We've harnessed the power of geometry to shape today's most advanced orthokeratology lens design for myopia management. We call it Euclid Emerald. It's orthokeratology elevated. The smart, modern way to offer Ortho-K.

Discover some of the unique advantages of our proprietary MyoK^{*} **lens design.**



Dual Alignment


gently align the peripheral cornea, offering better centration for a closer, comfortable fit

Empirical and Intuitive Design

offers control over the size of the treatment zone

MyoK ** Lens Design creates desirable peripheral myopic

creates desirable peripheral myopic defocus by ensuring plus power as close to the pupil as possible



Build and Differentiate Your Practice as a Leader in Proactive Myopia Management

Patients rely on your expertise and continuous care to achieve the best possible outcomes.

With today's growing myopia epidemic, Euclid Emerald Ortho-K treatment offers your practice the opportunity to be a part of

the exponential growth in myopia management for children.

Myopia is a Growing Opportunity

Ortho-K is an In-Demand Eyecare Specialty No Trial Lenses Required





Patient-Specific Lenses with Proprietary MyoK[™] Design

Each Euclid Emerald lens* is specifically designed for the individual patient's eyes. With Euclid Emerald, you don't select a lens that's the closest match from a fitting set. You provide patients the correct lens. It's tailored myopia management that complements your patient-centered practice.

Thanks to the power of three factor fitting, all we need are the following parameters: the Rx, the K readings and HVID. This efficient, empirical process is devoid of quesswork and yields category-leading

95.5% fit success:

Rx K's HVID

*Euclid Emerald Toric is available for patients with astigmatism. **With up to 1 lens exchange.



Build and Differentiate Your Practice as a Leader in Proactive Myopia Management

Patients rely on your expertise and continuous care to achieve the best possible outcomes.

For independent eye care professionals looking to differentiate **themselves, Euclid Emerald is the right product, in the right**

themselves, Euclid Emerald is the right product, in the right place, at the right time.

Euclid Emerald is Not Available Online

Easy Three Step Process to Get Started



Three Steps to Build Your Myopia Management **Practice with Euclid Emerald Ortho-K**

Get Emerald certified at getemeraldcertified.com

Schedule patient exam and record three parameters

Rx

K's

HVID

Send patient's parameters to sales@spcinternational.com or call (470) 208-7030

Our World-Class Consultation

Team is With You All the Way

Our professional consultation team makes orthokeratology easy for you. They're with you every step of the way to help you and your patients achieve the best possible outcomes. Euclid creates your patient-specific MyoK** lenses and our Professional Consultation Team helps make every step of the treatment process effortless.



The Global Leader in Myopia Management and Advanced Orthokeratology