# S P E C T R U M



DESIGNED BY **DATA**. DESIGNED TO **FIT**. **DESIGNED FOR LIFE** 

## S P E C T R U M



BostonSight SCLERAL is the first evidence-based lens design

that has leveraged patient fitting data and clinical experience - from 6 years

of patient outcomes on more than 7,000 eyes. BostonSight SCLERAL

helps you achieve a smart, predictable, and efficient fit with an optimal

#### ALIGNED FIT, RIGHT FROM THE START



90% of fits achieve optimal vision with the built-in front

surface eccentricity.



75% of fits achieved with the go-to Standard lens from the fitting set.




60% of fits don't require any change in haptic design.

Patients say it best - over 90% surveyed would recommend BostonSight SCLERAL to others.



Each fitting uses a single starting point and incorporates a right and left eye anatomical design. BostonSight SCLERAL allows you to customize each lens for every patient to ensure predictable and accurate outcomes.

### DISTINCT ADVANTAGES OF BOSTONSIGHT SCLERAL:



#### SplineCurve<sup>™</sup> теснлоlogy

Allows you to control the changes you want to make

#### **Quadrant-Specific Toric** PERIPHERAL HAPTIC SYSTEM (PHS™)

Allows you to control the changes you want to make in a quadrant-specific manner

#### SmartChannel<sup>™</sup> тесниоLOGY

Facilitating fits over anatomical obstacles and reducing suction

#### Front Surface ECCENTRICITY (FSE)

Designed to offer superior vision

#### Lens Diameters

16.0mm, 16.5mm, and 17.0mm

#### Comprehensive 28 Lens DIAGNOSTIC SET

Enhancing an efficient fitting process with right and left eye anatomical designs

Patients report a 60% increase in satisfaction when fit with BostonSight SCLERAL over prior lens.



### Comprehensive **FOR AN ENHANCED PRACTITIONER EXPERIENCE**

**Superior Value** The BostonSight SCLERAL fitting system offers superior value for practitioners – providing you with everything you need:

- No surprises or extra fees throughout the fitting process no extra charge for changes in diameter, base curve, material, toric periphery, front surface toricity, or channels.
- No special equipment necessary.



**FitAcademy™** With over 25 years of experience designing, manufacturing, contributing to peer-review literature, training clinicians and educating about scleral lenses, the BostonSight team is a leader and pioneer in the industry. The BostonSight SCLERAL team has developed FitAcademy – live and virtual workshops, designed to give clinicians an opportunity to improve scleral lens fitting, assessment, and management.

"I really appreciate the all-inclusive fee structure that BostonSight offers. I love not having to worry about surprise charges that may impact my patient. I also love the charitable side of BostonSight." - Dr. Karen Lee, Houston, TX



### S P E C T R U M



#### DESIGNED BY DATA. DESIGNED TO FIT. DESIGNED FOR LIFE





"The best set of diagnostic lenses on the market. Philosophically, the scleral lens fitting method is straightforward, promoting success in the end result".

Dr. Robert Davis, Chicago, Il

#### BostonSight SCLERAL LENS PARAMETER AVAILABILITY

DIAMETERS	16.0mm, 16.5mm, and 17.0mm
SPHERE POWER	-20.00 Diopters to +20.00 Diopters
SAGITTAL HEIGHT	2.0mm to 6.0mm in 0.1mm (100µm) steps
PERIPHERAL HAPTIC SYSTEM (PHS™)	Customizable
OPTIONS	Quadrant-specific toric PHS™ Front surface eccentricity Front surface toric Rx SmartChannels™
PREFERRED MATERIAL	Optimum Extra, Optimum Extreme, Optimum Infinite, Boston® EQII, Boston® XO2



BostonSight is a world-renowned leader in the treatment of patients with complex corneal disease and a pioneer in scleral lens design and technology. Founded in 1992 as a nonprofit eye healthcare organization, we launched BostonSight SCLERAL in 2017 to expand global access to scleral lens technology.



### S P E C T R U M

**\ +1 (470) 208-7030** ≥ sales@spctinternational.com www.spctinternational.com