

3-MONTHLY REPLACEMENT

SAPHIR

SILICONE HYDROGEL



Saphir provides comfortable, healthy contact lens wear to patients accustomed to a 3-monthly lens replacement.


-  **SPHERIC**
-  **MULTIFOCAL**
-  **TORIC**
-  **MULTIFOCAL TORIC**

PARAMETERS

Base curves (mm)	6.80 to 9.80 (0.30)
Diameters (mm)	13.00 to 16.00 (0.50)
Spheres (D)	±30.00 (0.25)
Cylinders	-0.75 to -8.00 (0.25)
Axes (°)	All (5°)
Additions	0.50 to 4.00 (0.50) CD/CN

MATERIAL





Type	Filcon 5B (60) [75%]
DK (ISO 9913-1-1998)	60
DK/t (-3.00 D)	50
Water Content	75%
Central Thickness (-3.00 D)	0.12
Cof	0.04
Modulus	0.29
Handling tint	Blue
Pack size	Single and 2-pack
Manufacturing Process	Lathed



Need fitting advice?

Check the Fitting Guide section of the catalogue.

OPTICAL DESIGN



Lens design parameters may change depending on the power