



High Definition Vision Aberration Control

 Silicone Hydrogel
  Monthly

 Spheric
  Toric



GrabSoo Technology

Patented technology for advanced and innovative function. Hydrophilic polymers grip tightly to each other creating a perfect contact lens surface, generating a clearer, stronger surface and helping to prevent protein deposits. In addition, blocks of highly lubricated material are the benefits of using **GrabSoo Technology**.

Parameters	Spheric	Toric
Material	Silicone Hydrogel	
Water Content	45%	
Base Curve	8.60 mm	
Diameter	14.20 mm	
Power	+6.00 to -6.00 (0.25 step) -6.50 to -12.00 (0.50 step) +6.50 to +8.00 (0.50 step)	SPH: +6.00 to -6.00 (0.25 step) -6.50 to -12.00 (0.50 step)
Cylinder Power	—	-0.75, -1.25, -1.75, -2.25, -2.75
Axis	—	10° to 180° (in 10° steps)
Dk/t (@-3.00D)	119	75
Dk	70	
Modulus	0.4 - 0.5	
Pack	6 lenses	

S P E C T R U M
INTERNATIONAL



✉ sales@spctinternational.com

☎ +1 (470) 208-7030

🌐 www.spctinternational.com

1345 Seaboard Industrial Blvd NW. Unit # C-2
Atlanta, Georgia. USA 