

SPECTRUM  
INTERNATIONAL

SPECTRUM  
color

# SPECTRUM color

**Multiple curves and high transparent materials helps the vision stay clear**

**High moisture, smooth tear circulation**


**Long-lasting comfort to wearers**

**SPECTRUM Color** brand contact lenses offer long-lasting comfort to wearers. Our lenses are manufactured with classic polymacon material, which is still widely recognized for its durable, dependable performance across all environments.


Advanced “in-toric” applied design reduces axis rotation. Circle “zone” gives superior visual correction. Toric lenses limitless optical range, cylinder and axis give wearers an individually, “PRECISE FIT”. The superior prism ballast design creates great visual acuity.

Securing visual clarity with  
a wide optical zone

The optical zone, formed by a sophisticated design, enables wearers to catch wide views from every angle.



Other Brand



SPECTRUM color

# SPECTRUM color

Material	62% polycarbon, 38% water content
Base Curve	8.60
Diameter	14.50mm
Graphic Color	Custom make
Power Range	<b>0.00 to ±25.00D</b> 0.00 to ±6.00D (0.25D steps) ±6.50 to ±25.00D (0.50D steps)
Cylinder Power	-0.50 to -7.00D (0.25D steps)
Axis	5° to 180° (in 5° steps)
Axis Marks	Three lines indication at 6 o'clock
Nominal Centre Thickness	0.130mm @-3.00D
Dk/t @-3.00D	7.51x10 <sup>-9</sup> (cm/s)(mLO <sub>2</sub> /ml.mmHg)
Oxygen Permeability Dk@35°C	9.77x10 <sup>-11</sup> (cm <sup>2</sup> /s)(mLO <sub>2</sub> /ml.mmHg)
Standard Pack Size	Vial: 1 lens


\*More colors available

\*\*We can design the color you require.




Yearly


## 3-tone color lens




**Tangerine**  
(CL365N)




**Olive**  
(CL370N)



**Turquoise**  
(CL304N)




**Gray**  
(CL301N)




**Dusk Blue**  
(CL367N)


Our high quality corrective lens gives you the  
clean and clear vision enhancing color lens



<b>Material</b>	Hioxifilcon A, 58% H <sub>2</sub> O
<b>Base Curve</b>	8.60
<b>Diameter</b>	14.20 mm
<b>Power Range</b>	<b>-10.00 to +6.00D</b> -0.50 to -6.00D (0.25D steps) -6.50 to -10.00D (0.50D steps) +0.50 to +4.00D (0.25D steps) +4.50 to +6.00D (0.50D steps)
<b>Packing Type</b>	10 Lenses
<b>Replacement</b>	Daily




**Light Blue**  
(CL240)




**Gray**  
(FL331)




**Hazel**  
(FL332)



**Blue**  
(FL333)




**Green**  
(FL334)




**Quarterly**


<b>Material</b>	Methafilcon A, 55% H <sub>2</sub> O
<b>Base Curve</b>	8.60
<b>Diameter</b>	14.20 mm
<b>Power Range</b>	<b>-10.00 to +6.00D</b> -0.50 to -6.00D (0.25D steps) -6.50 to -10.00D (0.50D steps) +0.25 to +4.00D (0.25D steps) +4.50 to +6.00D (0.50D steps)
<b>Packing Type</b>	2 Lenses
<b>Replacement</b>	Quarterly




**Tangerine**  
(FL301N)




**Violet**  
(FL302N)




**Emerald**  
(FL303N)




**Turquoise**  
(FL304N)




**Dusk Blue**  
(FL305N)




**Brown**  
(FL306N)




**Gray**  
(FL307N)




**Green**  
(FL308N)




**Marigold**  
(FL309N)





**Serenity**  
(FL310N)



**Cloud**  
(FL311N)



**Olive**  
(FL312N)



## Quarterly

<b>Material</b>	Polymacon, 38% H <sub>2</sub> O		
<b>Base Curve</b>	8.60	<b>Diameter</b>	14.20 mm
<b>Power Range</b>	<b>-10.00 to +6.00</b> -0.00 to -6.00 (0.25D steps) -6.50 to -10.00 (0.50D steps) +0.50 to +6.00 (0.50D steps)		
<b>Pack</b>	2 Lenses	<b>Replacement</b>	Quarterly



**Light Blue**  
(CL134N)




**Brown**  
(CL135N)



**Green**  
(CL136N)




**Violet**  
(CL137N)




**Aqua Blue**  
(CL138N)




**Gray**  
(CL139N)




**Gray**  
(CL301N)




**Brown**  
(CL302N)




**Violet**  
(CL303N)




**Turquoise**  
(CL304N)



**Green**  
(CL305N)




**Marigold**  
(CL364N)




**Tangerine**  
(CL365N)




**Emerald**  
(CL366N)




**Dusk Blue**  
(CL367N)



**Serenity**  
(CL368N)



**Cloud**  
(CL369N)



**Olive**  
(CL370N)

*\*More colors available*

*\*\*We can design the color you require.*



## Spectrum Color Halloween 1 Day

<b>Material</b>	Polymacon, 38% H <sub>2</sub> O		
<b>Base Curve</b>	8.60	<b>Diameter</b>	14.20 mm
<b>Power Range</b>	0.00 to -6.00D (0.25D steps)		
<b>Packing Type</b>	2 Lenses	<b>Replacement</b>	Daily



**Black Out**  
(FN101N)




**Red Vampire**  
(FN102N)




**White Out**  
(FN103N)



**Angelic White**  
(FN104N)



**Blue Walker**  
(FN105N)



**Zombie Yellow**  
(FN106N)

## Spectrum Color Halloween

<b>Material</b>	Polymacon, 38% H <sub>2</sub> O		
<b>Base Curve</b>	8.60	<b>Diameter</b>	14.20 mm
<b>Power Range</b>	0.00		
<b>Packing Type</b>	2 Lenses	<b>Replacement</b>	Monthly



**Angelic Blue**  
(FN002N)



**Angelic Red**  
(FN004N)



**Black Out**  
(FN009N)



**Black Wolf**  
(FN010)




**Lestat**  
(FN030N)



**White Out**  
(FN067N)



**Banshee**  
(FN006)



**Manson**  
(FN032)



**Red Cat**  
(FN042N)



**Red Vampire**  
(FN044)

# SPECTRUM

INTERNATIONAL



✉ [sales@spctinternational.com](mailto:sales@spctinternational.com)

☎ +1 (470) 208-7030



[www.spctinternational.com](http://www.spctinternational.com)