## Multipurpose Solutions





### Oté Sensation All-in-one

Premium hyaluronate solution for soft contact lenses

In addition to complete protein removal, this innovative hyaluronic solution has a dual antibacterial action. With the addition of high-quality hyaluronate, Oté Sensation offers incredibly high wearing comfort.

#### What are its functions?

Clean, insert, condition, hydrate and remove the protein layer, disinfect, rinse and moisten all types of soft lenses.

## **Oté Fine**

# Multipurpose solution all in one

Comfort solution for hard gas-permeable contact lenses.


This all-in-one solution for RGP contact lenses provides optimal cleaning and comfort. This fluid also has very good disinfecting qualities and is ideal for Ortho-k contact lenses.

#### What are its functions?

Cleaning, disinfection, conservation, placement and wetting of RGP contact lenses.

#### SPECTRUM

INTERNATIONAL



✓ sales@spctinternational.com



